

Product Datasheet

Recombinant Mouse IL-23R/IL-23 Receptor Protein, Fc His Tag (orb1550707)

Description	Recombinant Mouse IL-23R/IL-23 Receptor Protein, Fc His Tag
Species/Host	Human
Reactivity	Mouse
Endotoxins	< 1 EU/ug by LAL method.
Target	IL23R
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1782, Recombinant Mouse IL-23R/IL-23 Receptor Protein is produced by Human Cells expression system. The target protein is expressed with sequence (Gly24-Asp372) of mouse IL-23R/IL-23 Receptor (Accession #Q8K4B4) fused with an Fc tag at the C-terminus.
Expiration Date	6 months from date of receipt.