

Product Datasheet

Recombinant Human IL-18RAcP Protein, His Tag (orb1550715)

Description	Recombinant Human IL-18RAcP Protein, His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 1.0 EU/ug by LAL method.
Target	IL18RAcP
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1774, Recombinant Human IL-18RAcP Protein is produced by Mammalian expression system. The target protein is expressed with sequence (Phe20-Arg356) of human IL-18RAcP (Accession #O95256) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.