

Product Datasheet

Recombinant Human CD40L Protein (orb1550716)

Description	Recombinant Human CD40L Protein
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 1.0 EU/ug by LAL method.
Target	CD40L
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1773, Recombinant Human CD40L Protein is produced by E. coli expression system. The target protein is expressed with sequence (Met113-Leu261) of human CD40L (Accession #P29965) fused with no tags.
Expiration Date	6 months from date of receipt.