

Product Datasheet

Recombinant Human LILRB1 Protein, His Tag (orb1550720)

Description	Recombinant Human LILRB1 Protein, His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 1.0 EU/ug by LAL method.
Target	LILRB1
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1769, Recombinant Human LILRB1 Protein is produced by Mammalian expression system. The target protein is expressed with sequence (Gly24-His Tag458) of human LILRB1 (Accession #Q8NHL6) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.