

## Product Datasheet

# Recombinant Human FCGRT&B2M Protein, His Tag (orb1550721)

<b>Description</b>	Recombinant Human FCGRT&B2M Protein, His Tag
<b>Species/Host</b>	Human
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 0.1 EU/ug by LAL method.
<b>Target</b>	FCGRT
<b>Note</b>	For research use only
<b>Purity</b>	> 97% by SDS-PAGE.
<b>Source</b>	KMP1768, Recombinant Human FCGRT&B2M Protein is produced by HEK293 Cells expression system. The target heterodimer proteins are co-expressed with the sequence (Ala24-Ser297) of human FCGRT (Accession #NP_004098.1) fused with a 6×His Tag at the C-terminus and the sequence (Ile21-Met119) of human B2M (Accession #NP_004039.1).
<b>Expiration Date</b>	6 months from date of receipt.