

## Product Datasheet

# Recombinant Mouse PLGF/PGF Protein, Fc & His Tag (orb1550725)

|                        |  |
|------------------------|--|
| <b>Description</b>     | Recombinant Mouse PLGF/PGF Protein, Fc & His Tag   |
| <b>Species/Host</b>    | Human  |
| <b>Reactivity</b>      | Mouse  |
| <b>Endotoxins</b>      | < 1.0 EU/ug by LAL method.   |
| <b>Target</b>          | PLGF   |
| <b>Note</b>            | For research use only  |
| <b>Purity</b>          | > 87% by SDS-PAGE.   |
| <b>Source</b>          | KMP1764, Recombinant Mouse PLGF/PGF Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Met1-Pro158) of mouse PLGF/PGF (Accession #NP_032853.1) fused with an Fc, 6×His Tag at the C-terminus. |
| <b>Expiration Date</b> | 6 months from date of receipt.   |