

Product Datasheet

Recombinant Human PMP2 Protein, His Tag (orb1550726)

Description	Recombinant Human PMP2 Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 1.0 EU/ug by LAL method.
Target	PMP2
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1763, Recombinant Human PMP2 Protein is produced by E. coli expression system. The target protein is expressed with sequence (Met1-Val132) of human PMP2 (Accession #P02689) fused with a 6×His Tag at the N-terminus.
Expiration Date	6 months from date of receipt.