

## Product Datasheet

### Recombinant Human IL-1F10 Protein (orb1550728)

<b>Description</b>	Recombinant Human IL-1F10 Protein
<b>Species/Host</b>	E. coli
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1.0 EU/ug by LAL method.
<b>Target</b>	IL1F10
<b>Note</b>	For research use only
<b>Purity</b>	> 95% by SDS-PAGE.
<b>Source</b>	KMP1761, Recombinant Human IL-1F10 Protein is produced by E. coli expression system. The target protein is expressed with sequence (Met1-Trp152) of human IL-1F10 (Accession #Q8WWZ1) fused with no tags.
<b>Expiration Date</b>	6 months from date of receipt.