

Product Datasheet

Recombinant Mouse IL-21R/Interleukin-21 Receptor Protein, His Tag (orb1550730)

Description	Recombinant Mouse IL-21R/Interleukin-21 Receptor Protein, His Tag
Species/Host	Human
Reactivity	Mouse
Endotoxins	< 0.1 EU/ug by LAL method.
Target	IL21R
Note	For research use only
Purity	> 97% by SDS-PAGE.
Source	KMP1759, Recombinant Mouse IL-21R/Interleukin-21 Receptor Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Cys20-Pro236) of mouse IL-21R/Interleukin-21 Receptor (Accession #NP_068687.1) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.