

## Product Datasheet

# Recombinant Human Adiponectin/Acrp30/ADIPOQ Protein, His Tag (orb1550731)

<b>Description</b>	Recombinant Human Adiponectin/Acrp30/ADIPOQ Protein, His Tag
<b>Species/Host</b>	Human
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1 EU/ug by LAL method.
<b>Target</b>	Adiponectin
<b>Note</b>	For research use only
<b>Purity</b>	> 95% by SDS-PAGE.
<b>Source</b>	KMP1758, Recombinant Human Adiponectin/Acrp30/ADIPOQ Protein is produced by Human Cells expression system. The target protein is expressed with sequence (Glu19-Asn244) of human Adiponectin/Acrp30/ADIPOQ (Accession #Q15848) fused with a 6xHis Tag at the C-terminus.
<b>Expiration Date</b>	6 months from date of receipt.