

Product Datasheet

Recombinant Mouse IL-7 Protein, His Tag (orb1550732)

Description	Recombinant Mouse IL-7 Protein, His Tag
Species/Host	Mammal
Reactivity	Mouse
Endotoxins	< 1.0 EU/ug by LAL method.
Target	IL7
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1757, Recombinant Mouse IL-7 Protein is produced by Mammalian expression system. The target protein is expressed with sequence (Glu26-Ile154) of mouse IL-7 (Accession #P10168) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.