

Product Datasheet

Recombinant Human IL-23R Protein, His Tag (orb1550739)

Description	Recombinant Human IL-23R Protein, His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 1.0 EU/ug by LAL method.
Target	IL23R
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1750, Recombinant Human IL-23R Protein is produced by Mammalian expression system. The target protein is expressed with sequence (Gly24-Asp353) of human IL-23R (Accession #Q5VWK5) fused with a 6xHis Tag at the C-terminus.
Expiration Date	6 months from date of receipt.