

Product Datasheet

Recombinant Human Vitronectin/VTN/VN Protein, His Tag (orb1550746)

Description	Recombinant Human Vitronectin/VTN/VN Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	Vitronectin
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1743, Recombinant Human Vitronectin/VTN/VN Protein is produced by Human Cells expression system. The target protein is expressed with sequence (Asp20-Leu478) of human Vitronectin/VTN/VN (Accession #P04004) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.