

Product Datasheet

Recombinant Human PDGF-BB Protein (orb1550747)

Description	Recombinant Human PDGF-BB Protein
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	PDGFBB
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1742, Recombinant Human PDGF-BB Protein is produced by E. coli expression system. The target protein is expressed with sequence (Ser82-Thr190) of human PDGF-BB (Accession #P01127) fused with an initial Met at the N-terminus.
Expiration Date	6 months from date of receipt.