

Product Datasheet

Recombinant Human CXCL7/NAP-2 Protein (orb1550754)

Description	Recombinant Human CXCL7/NAP-2 Protein
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	CXCL7
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1735, Recombinant Human CXCL7/NAP-2 Protein is produced by E. coli expression system. The target protein is expressed with sequence (Ala59-Asp128) of human CXCL7/NAP-2 (Accession #P02775).
Expiration Date	6 months from date of receipt.