

## Product Datasheet

### Recombinant Human BMP-2 Protein (orb1550756)

<b>Description</b>	Recombinant Human BMP-2 Protein
<b>Species/Host</b>	E. coli
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1 EU/ug by LAL method.
<b>Target</b>	BMP2
<b>Note</b>	For research use only
<b>Purity</b>	> 95% by SDS-PAGE.
<b>Source</b>	KMP1733, Recombinant Human BMP-2 Protein is produced by E. coli expression system. The target protein is expressed with sequence (Gln283-Arg396) of human BMP-2 (Accession #P12643) fused with an initial Met at the N-terminus.
<b>Expiration Date</b>	6 months from date of receipt.