

## Product Datasheet

### Recombinant Human IgG4 Fc Protein (orb1550757)

<b>Description</b>	Recombinant Human IgG4 Fc Protein
<b>Species/Host</b>	Human
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1 EU/ug by LAL method.
<b>Target</b>	IgG4 Fc
<b>Note</b>	For research use only
<b>Purity</b>	> 95% by SDS-PAGE.
<b>Source</b>	KMP1732, Recombinant Human IgG4 Fc Protein is produced by Human Cells expression system. The target protein is expressed with sequence (Glu99-Gly326) of human IgG4 Fc (Accession #P01861).
<b>Expiration Date</b>	6 months from date of receipt.