

Product Datasheet

Recombinant Mouse B7-H3/CD276 Protein, His Tag (orb1550759)

Description	Recombinant Mouse B7-H3/CD276 Protein, His Tag
Species/Host	Human
Reactivity	Mouse
Endotoxins	< 1 EU/ug by LAL method.
Target	B7H3
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1730, Recombinant Mouse B7-H3/CD276 Protein is produced by Human Cells expression system. The target protein is expressed with sequence (Val29-Phe244) of mouse B7-H3/CD276 (Accession #Q8VE98) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.