

Product Datasheet

Recombinant Mouse TNFSF4/OX40 Ligand Protein, His Tag (orb1550761)

Description	Recombinant Mouse TNFSF4/OX40 Ligand Protein, His Tag
Species/Host	Human
Reactivity	Mouse
Endotoxins	< 1 EU/ug by LAL method.
Target	TNFSF4
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1728, Recombinant Mouse TNFSF4/OX40 Ligand Protein is produced by Human Cells expression system. The target protein is expressed with sequence (Ser51-Leu198) of mouse TNFSF4/OX40 Ligand (Accession #P43488) fused with an 8xHis Tag at the N-terminus.
Expiration Date	6 months from date of receipt.