

Product Datasheet

Recombinant Human CD94/Natural Killer Cells Antigen CD94 Protein, His Tag (orb1550766)

Description	Recombinant Human CD94/Natural Killer Cells Antigen CD94 Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	CD94
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1723, Recombinant Human CD94/Natural Killer Cells Antigen CD94 Protein is produced by Human Cells expression system. The target protein is expressed with sequence (Ser34-Ile179) of human CD94/Natural Killer Cells Antigen CD94 (Accession #Q13241) fused with a 6×His Tag at the N-terminus.
Expiration Date	6 months from date of receipt.