

Product Datasheet

Recombinant Human CD160 Protein, His Tag (orb1550767)

Description	Recombinant Human CD160 Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	CD160
Note	For research use only
Purity	> 90% by SDS-PAGE.
Source	KMP1722, Recombinant Human CD160 Protein is produced by Human Cells expression system. The target protein is expressed with sequence (Ile27-Ser159) of human CD160 (Accession #O95971) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.