

Product Datasheet

Recombinant Human CD79B/B29 Protein, His Tag (orb1550768)

Description	Recombinant Human CD79B/B29 Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	CD79B
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1721, Recombinant Human CD79B/B29 Protein is produced by Human Cells expression system. The target protein is expressed with sequence (Ala29-Asp159) of human CD79B/B29 (Accession #P40259) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.