

Product Datasheet

Recombinant Human Heat Shock 47kDa/SERPINH1 Protein, His Tag (orb1550770)

Description	Recombinant Human Heat Shock 47kDa/SERPINH1 Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	Heat Shock 47kDa
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1719, Recombinant Human Heat Shock 47kDa/SERPINH1 Protein is produced by Human Cells expression system. The target protein is expressed with sequence (Ala19-Leu418) of human Heat Shock 47kDa/SERPINH1 (Accession #P50454) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.