

Product Datasheet

Recombinant Mouse IL-10 Protein (orb1550772)

Description	Recombinant Mouse IL-10 Protein
Species/Host	E. coli
Reactivity	Mouse
Endotoxins	< 1 EU/ug by LAL method.
Target	IL10
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1717, Recombinant Mouse IL-10 Protein is produced by E. coli expression system. The target protein is expressed with sequence (Ser19-Ser178) of mouse IL-10 (Accession #P18893) fused with an initial Met at the N-terminus.
Expiration Date	6 months from date of receipt.