

Product Datasheet

Recombinant Human R-spondin-3/RSPO3 Protein, Fc & His Tag (orb1550774)

Description	Recombinant Human R-spondin-3/RSPO3 Protein, Fc & His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	R-spondin3
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1715, Recombinant Human R-spondin-3/RSPO3 Protein is produced by Human Cell expression system. The target protein is expressed with sequence (Gln22-His Tag272) of human R-spondin-3/RSPO3 (Accession #Q9BXY4) fused with an Fc, 6xHis Tag at the C-terminus.
Expiration Date	6 months from date of receipt.