

Product Datasheet

Recombinant Human S100A8/Protein S100-A8 Protein, His Tag (orb1550775)

Description	Recombinant Human S100A8/Protein S100-A8 Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	S100A8
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1714, Recombinant Human S100A8/Protein S100-A8 Protein is produced by E. coli expression system. The target protein is expressed with sequence (Met1-Glu93) of human S100A8/Protein S100-A8 (Accession #P05109) fused with a 6xHis Tag at the C-terminus.
Expiration Date	6 months from date of receipt.