

Product Datasheet

Recombinant Human Prolactin/PRL Protein, His Tag (orb1550776)

Description	Recombinant Human Prolactin/PRL Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	Prolactin
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1713, Recombinant Human Prolactin/PRL Protein is produced by Human Cell expression system. The target protein is expressed with sequence (Leu29-Cys227) of human Prolactin/PRL (Accession #P01236) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.