

Product Datasheet

Recombinant Mouse Pro-EGF Protein, His Tag (orb1550779)

Description	Recombinant Mouse Pro-EGF Protein, His Tag
Species/Host	Human
Reactivity	Mouse
Endotoxins	< 0.1 EU/ug by LAL method.
Target	ProEGF
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1710, Recombinant Mouse Pro-EGF Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Trp29-Arg1029) of mouse Pro-EGF (Accession #NP_034243.2) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.