

Product Datasheet

Recombinant Mouse TNFR1/CD120a/TNFRSF1A Protein, His Tag (orb1550785)

Description	Recombinant Mouse TNFR1/CD120a/TNFRSF1A Protein, His Tag
Species/Host	Human
Reactivity	Mouse
Endotoxins	< 0.1 EU/ug by LAL method.
Target	TNFR1
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1704, Recombinant Mouse TNFR1/CD120a/TNFRSF1A Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Met1-Ala212) of mouse TNFR1/CD120a/TNFRSF1A (Accession #NP_035739.2) fused with a 6xHis Tag at the C-terminus.
Expiration Date	6 months from date of receipt.