

Product Datasheet

Recombinant Mouse CXCL16 Protein, His Tag (orb1550792)

Description	Recombinant Mouse CXCL16 Protein, His Tag
Species/Host	Human
Reactivity	Mouse
Endotoxins	< 0.1 EU/ug by LAL method.
Target	CXCL16
Note	For research use only
Purity	> 97% by SDS-PAGE.
Source	KMP1697, Recombinant Mouse CXCL16 Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Asn27-Trp201) of mouse CXCL16 (Accession #NP_075647.3) fused with a 6×His Tag at the C- terminus.
Expiration Date	6 months from date of receipt.