

Product Datasheet

Recombinant human Nectin-1/PVRL1/CD111 Protein, Fc & His Tag (orb1550794)

Description	Recombinant human Nectin-1/PVRL1/CD111 Protein, Fc & His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	Nectin1
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1695, Recombinant human Nectin-1/PVRL1/CD111 Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Met1-Thr334) of human Nectin-1/PVRL1/CD111 (Accession #NP_002846.3) fused with an Fc, 6xHis Tag at the C-terminus.
Expiration Date	6 months from date of receipt.