

Product Datasheet

Recombinant Mouse Dkk-1 Protein, His Tag (orb1550796)

Description	Recombinant Mouse Dkk-1 Protein, His Tag
Species/Host	Human
Reactivity	Mouse
Endotoxins	< 0.1 EU/ug by LAL method.
Target	Dkk1
Note	For research use only
Purity	> 87% by SDS-PAGE.
Source	KMP1693, Recombinant Mouse Dkk-1 Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Met1-His Tag272) of mouse Dkk-1 (Accession #NP_034181.2.) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.