

Product Datasheet

Recombinant Human TWSG1 Protein, His Tag (orb1550804)

Description	Recombinant Human TWSG1 Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	TWSG1
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1685, Recombinant Human TWSG1 Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Met1-Phe223) of human TWSG1 (Accession #NP_065699.1) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.