

Product Datasheet

Recombinant Human PDGFAA/PDGF-AA Protein (orb1550807)

Description	Recombinant Human PDGFAA/PDGF-AA Protein
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	PDGFAA
Note	For research use only
Purity	> 90% by SDS-PAGE.
Source	KMP1682, Recombinant Human PDGFAA/PDGF-AA Protein is produced by E. coli expression system. The target protein is expressed with sequence (Ser87-Thr211) of human PDGFAA/PDGF-AA (Accession #P04085).
Expiration Date	6 months from date of receipt.