

Product Datasheet

Recombinant Human IL-22 Protein (orb1550808)

Description	Recombinant Human IL-22 Protein
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	IL22
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1681, Recombinant Human IL-22 Protein is produced by E. coli expression system. The target protein is expressed with sequence (Ala34-Ile179) of human IL-22 (Accession #Q9GZX6) fused with an initial Met at the N-terminus.
Expiration Date	6 months from date of receipt.