

Product Datasheet

Recombinant Human Jagged-1/JAG1 Protein, Fc His Tag (orb1550812)

Description	Recombinant Human Jagged-1/JAG1 Protein, Fc His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	Jagged1
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1677, Recombinant Human Jagged-1/JAG1 Protein is produced by Human Cells expression system. The target protein is expressed with sequence (Gln34-Ser1046) of human Jagged-1/JAG1 (Accession #P78504) fused with an Fc tag at the C-terminus.
Expiration Date	6 months from date of receipt.