

Product Datasheet

Recombinant Human Fc epsilon RII/CD23 Protein, His Tag (orb1550820)

Description	Recombinant Human Fc epsilon RII/CD23 Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	Fc epsilon RII
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1669, Recombinant Human Fc epsilon RII/CD23 Protein is produced by Human Cells expression system. The target protein is expressed with sequence (Asp48-Ser321) of human Fc epsilon RII/CD23 (Accession #P06734) fused with an 8×His Tag at the N-terminus.
Expiration Date	6 months from date of receipt.