

Product Datasheet

Recombinant Mouse TNFSF11/RANK L/TRANCE Protein (orb1550823)

Description	Recombinant Mouse TNFSF11/RANK L/TRANCE Protein
Species/Host	E. coli
Reactivity	Mouse
Endotoxins	< 1 EU/ug by LAL method.
Target	TNFSF11
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1666, Recombinant Mouse TNFSF11/RANK L/TRANCE Protein is produced by E. coli expression system. The target protein is expressed with sequence (Gln137-Asp316) of mouse TNFSF11/RANK L/TRANCE (Accession #O35235).
Expiration Date	6 months from date of receipt.