

Product Datasheet

Recombinant Mouse CD155/PVR Protein, His Tag (orb1550824)

Description	Recombinant Mouse CD155/PVR Protein, His Tag
Species/Host	Human
Reactivity	Mouse
Endotoxins	< 1 EU/ug by LAL method.
Target	CD155
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1665, Recombinant Mouse CD155/PVR Protein is produced by Human Cells expression system. The target protein is expressed with sequence (Asp29-Leu348) of mouse CD155/PVR (Accession #Q8K094) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.