

Product Datasheet

Recombinant Human IL-13RA1/IL-13 R alpha 1 Protein, His Tag (orb1550825)

Description	Recombinant Human IL-13RA1/IL-13 R alpha 1 Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	IL13RA1
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1664, Recombinant Human IL-13RA1/IL-13 R alpha 1 Protein is produced by Human Cells expression system. The target protein is expressed with sequence (Gly22-Thr343) of human IL-13RA1/IL-13 R alpha 1 (Accession #P78552) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.