

Product Datasheet

Recombinant Mouse Complement Component C3a/C3a Protein (orb1550828)

Description	Recombinant Mouse Complement Component C3a/C3a Protein
Species/Host	E. coli
Reactivity	Mouse
Endotoxins	< 1 EU/ug by LAL method.
Target	Complement Component C3a
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1661, Recombinant Mouse Complement Component C3a/C3a Protein is produced by E. coli expression system. The target protein is expressed with sequence (Ser671-Arg748) of mouse Complement Component C3a/C3a (Accession #P01027).
Expiration Date	6 months from date of receipt.