

Product Datasheet

Recombinant Mouse IL-18/IL-1F4 Protein, His Tag (orb1550836)

Description	Recombinant Mouse IL-18/IL-1F4 Protein, His Tag
Reactivity	Mouse
Endotoxins	> 95% by SDS-PAGE.
Target	IL18
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	Recombinant Mouse IL-18/IL-1F4 Protein is produced by E. coli expression system. The target protein is expressed with sequence (Asn36-Ser192) of mouse IL-18/IL-1F4 (Accession #P70380) fused with an initial Met at the N-terminus and a 6His Tag at the N-terminus.
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)