

Product Datasheet

Recombinant Human TGF-beta 3/TGFB3 Protein (orb1550837)

Description	Recombinant Human TGF-beta 3/TGFB3 Protein
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	TGFBeta 3
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1652, Recombinant Human TGF-beta 3/TGFB3 Protein is produced by Human Cells expression system. The target protein is expressed with sequence (Ala301-Ser412(Tyr340Phe)) of human TGF-beta 3/TGFB3 (Accession #P10600).
Expiration Date	6 months from date of receipt.