

Product Datasheet

Recombinant Mouse Shh/Sonic Hedgehog Protein (orb1550841)

Description	Recombinant Mouse Shh/Sonic Hedgehog Protein
Species/Host	E. coli
Reactivity	Mouse
Endotoxins	< 1 EU/ug by LAL method.
Target	Shh
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1648, Recombinant Mouse Shh/Sonic Hedgehog Protein is produced by E. coli expression system. The target protein is expressed with sequence (Cys25-Gly198) of mouse Shh/Sonic Hedgehog (Accession #Q62226) fused with an initial Met at the N-terminus.
Expiration Date	6 months from date of receipt.