

## Product Datasheet

# Recombinant Human IL-20RA/IL-20 R alpha Protein, His Tag (orb1550850)

<b>Description</b>	Recombinant Human IL-20RA/IL-20 R alpha Protein, His Tag
<b>Species/Host</b>	Human
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1 EU/ug by LAL method.
<b>Target</b>	IL20RA
<b>Note</b>	For research use only
<b>Purity</b>	> 95% by SDS-PAGE.
<b>Source</b>	KMP1639, Recombinant Human IL-20RA/IL-20 R alpha Protein is produced by Human Cells expression system. The target protein is expressed with sequence (Val30-Lys250) of human IL-20RA/IL-20 R alpha (Accession #Q9UHF4) fused with a 6×His Tag at the C-terminus.
<b>Expiration Date</b>	6 months from date of receipt.