

Product Datasheet

Recombinant Human CD157/BST1/ADP-ribosyl cyclase 2/Cyclic ADP-ribose hydrolase 2 Protein, His Tag (orb1550852)

Description	Recombinant Human CD157/BST1/ADP-ribosyl cyclase 2/Cyclic ADP-ribose hydrolase 2 Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	CD157
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1637, Recombinant Human CD157/BST1/ADP-ribosyl cyclase 2/Cyclic ADP-ribose hydrolase 2 Protein is produced by Human Cells expression system. The target protein is expressed with sequence (Gly29-Lys292) of human CD157/BST1/ADP-ribosyl cyclase 2/Cyclic ADP-ribose hydrolase 2 (Accession #Q10588) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.