

## Product Datasheet

# Recombinant Human VEGF-D Protein, His Tag (orb1550853)

<b>Description</b>	Recombinant Human VEGF-D Protein, His Tag
<b>Species/Host</b>	Human
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1 EU/ug by LAL method.
<b>Target</b>	VEGFD
<b>Note</b>	For research use only
<b>Purity</b>	> 95% by SDS-PAGE.
<b>Source</b>	KMP1636, Recombinant Human VEGF-D Protein is produced by Human Cells expression system. The target protein is expressed with sequence (Phe93-Ser201) of human VEGF-D (Accession #O43915) fused with a 6xHis Tag at the C-terminus.
<b>Expiration Date</b>	6 months from date of receipt.