

Product Datasheet

Recombinant Human Chitinase-3-like protein 1/CHI3L1 Protein, His Tag (orb1550854)

Description	Recombinant Human Chitinase-3-like protein 1/CHI3L1 Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	Chitinase3like 1
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1635, Recombinant Human Chitinase-3-like protein 1/CHI3L1 Protein is produced by Human Cells expression system. The target protein is expressed with sequence (Tyr22-Thr383) of human Chitinase-3-like protein 1/CHI3L1 (Accession #P36222) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.