

Product Datasheet

Recombinant Human IFN-alpha/beta R1 Protein, His Tag (orb1550855)

Description	Recombinant Human IFN-alpha/beta R1 Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	IFNalpha
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1634, Recombinant Human IFN-alpha/beta R1 Protein is produced by Human Cells expression system. The target protein is expressed with sequence (Lys28-Lys436) of human IFN-alpha/beta R1 (Accession #P17181) fused with a 6xHis Tag at the C-terminus.
Expiration Date	6 months from date of receipt.