

## Product Datasheet

# Recombinant Human NT-3/Neurotrophin-3 Protein (orb1550857)

<b>Description</b>	Recombinant Human NT-3/Neurotrophin-3 Protein
<b>Species/Host</b>	E. coli
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1 EU/ug by LAL method.
<b>Target</b>	NT3
<b>Note</b>	For research use only
<b>Purity</b>	> 95% by SDS-PAGE.
<b>Source</b>	KMP1632, Recombinant Human NT-3/Neurotrophin-3 Protein is produced by E. coli expression system. The target protein is expressed with sequence (Tyr139-Thr257) of human NT-3/Neurotrophin-3 (Accession #P20783).
<b>Expiration Date</b>	6 months from date of receipt.